



Water Distiller MLW-D10

Water Distiller MLW-D10

Features:

- * Stainless steel material, adopt stamping and welding technology.
- * Water shortage and power-off alarm; automatically make up water and heat again.
- * Power, heating indicator and show water shortage.

Technical Parameters:

Model no.	MLW-D10
Water Output	≥10L/H
Consumption	7.5KW
Power Supply	380V/50Hz
External Size(Dia.*H)	235*855mm
Package Size(W*D*H)	360*360*920mm
Gross Weight	9.5kg