



# ***Automatic Electronic Analytical Balance MLA-B15***

# Automatic Electronic Analytical Balance MLA-B15

## Products Description:

- \* Multiple functions guarantee easy and reliable operation, including auto-counting unit conversion (g, oz, ct...), full-scale faring, zero-memory etc
- \* Built-in RS232 automatic calibration guarantee direct connection to external equipment such as computer and printer
- \* Built-in automatic calibration guarantee easier, faster, stable and accurate

## Features:

- \* Humanized design, novel design and four-side transparent windshield
- \* High precision magnet sensor keeps the weighing more fast and stable
- \* Ultra-large backlight LCD display brings clearer and more comfortable visual sense to operate

## Technical Parameters:

|                   |                |
|-------------------|----------------|
| <b>Model no.</b>  | <b>MLA-B15</b> |
| Capacity          | 0~200g         |
| Readability       | 0.1mg          |
| Scale Size        | Φ80mm          |
| External Size     | 350*215*340mm  |
| Package Size      | 480*375*515mm  |
| G.W.              | 11KG           |
| Work Space Height | 230mm          |
| Calibration       | Internal Cal   |