



FORCED AIR DRYING OVEN MLL-H14

Forced Air Drying Oven MLL-H14

Features:

- * External material: cold-rolled steel with anti-bacteria powder coating.
- * Stainless steel inner chamber, round angle structure, adjustable shelves.
- * Double-layer glass observation window, adjustable airtight buckle lock.
- * PID microprocessor intelligent temperature controller.
- * Over-temperature protection.
- * Optional independent temperature protector

Technical Parameters:

Model no.	MLL-H14
Capacity	140L
Temp. Range	RT+10~250 /300°C
Temp. Precision	0.1°C
Temp. Fluctuation	±0.5°C
Ambient Temp.	5~40°C
Timing Range	1~9999min
Shelves No.	2 pcs
Display	LCD
Consumption	1500W
Power Supply	AC110/220V±10%,50/60Hz;
Internal Size (W*D*H)mm	450*557*520
External Size (W*D*H)mm	600*740*920
Packing Size (W*D*H)mm	730*860*1080
Gross Weight(kg)	95