



ULTRASONIC BATH

MUW-A40

Ultrasonic Bath MUW-A40

Introduction:

Ultrasonic cleaning is a process that uses ultrasonic and an appropriate cleaning solvent to clean items.

Features:

- * Standard basket and lid.
- * Fully microprocessor controller and programmable.
- * Digital display, easy operation for time and temperature setting.
- * Count-down time and actual temperature real time displaying.
- * Power Supply: AC220V±10%, 50/60Hz; 110V±10%, 60Hz

Technical Parameters:

Model No.	MUW-A40
Capacity (L)	10
Ultrasonic Frequency (Khz)	40
Time Range (min)	1~99
HeatingTemp. (°C)	20~80°C
Ultrasonic Power (W)	240
Heating Power (W)	250
Tank Size: L/W/H (mm)	300*240*150
External Size: L*W*H(mm)	330*270*270